


Product Name: Recombinant Human ATXN8OS Protein
Catalog #: PHH8365



Summary

Name	AT8OS_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P0DMR3
Host	HEK293
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FTSKSRPEP DMLYAQKAKG WQLTQDSGGW EVQDQCTRIW SKENLLALNT HSRRQKGKRE NKVCVSTWQK

image



Background

Alternative Names	Putative protein ATXN8OS (ATXN8 opposite strand) (Spinocerebellar ataxia 8) (kelch-like 1 antisense)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use