



Product Name: Recombinant Human PEX6 Protein
Catalog #: PHH8386



Summary

Name	PEX6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q13608
Host	HEK293
Species	Human
Molecular Mass	51 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KNSGLAGGLT EEDEGELCAA GFPLLAEDFG QALEQLQTAH SQAVGAPKIP SVSWHDVGGL QEVKKEILET IQLPLEHPEL

image



Background

Alternative Names	Peroxisome assembly factor 2 (PAF-2) (Peroxin-6) (Peroxisomal biogenesis factor 6) (Peroxisomal-type ATPase 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use