



Product Name: Recombinant Human CCDC6 Protein
Catalog #: PHH8396



Summary

Name	CCDC6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q16204
Host	HEK293
Species	Human
Molecular Mass	55 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GGGGKSGGIV ISPFRLLELT NRLASLQQEN KVLKIELETY KLKCKALQEE NRDLRKASVT IQARAEQEEE FISNTLFKKI QALQKEKETL AVNYEKEEEF LTNELSRKLM

image



Background

Alternative Names	Coiled-coil domain-containing protein 6 (Papillary thyroid carcinoma-encoded protein) (Protein H4)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use