



Product Name: Recombinant Human TAF15 Protein
Catalog #: PHH8414



Summary

Name	RBP56_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q92804
Host	HEK293
Species	Human
Molecular Mass	52 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LGEGVSTDQV GEFFKQIGII KTNKKTGKPM INLYTDKDTG KPKGATVSF DDPPSAKAAI DWFDGKEFHG NIIKVSFATR

image



Background

Alternative Names TATA-binding protein-associated factor 2N (68 kDa TATA-binding protein-associated factor) (TAF(II)68) (TAFII68) (RNA-binding protein 56)

Note

For Research Use Only , Not for Diagnostic Use