


Summary

Name	ELAV4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P26378
Host	HEK293
Species	Human
Molecular Mass	84 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MEWNLKMI STMEPQVSN PTSNTSNGPS SNNRNCPSM QTGATTDDSK TNLIVNYLPQ NMTQEEFRSL FGSIGEIESC KLVRDKITGQ SLGYGFVNYI DPKDAEKAIN TLNGLRLQTK TIKVSYARPS SASIRDANLY VSGLPKTMTQ KELEQLFSQY GRIITSRILV DQVTGVSRGV GFIRFDKRIE AEEAIKGLNG QKPSGATEPI TVKFANNPSQ KSSQALLSQL YQSPNRRYPG PLHHQAQFRF LDNLLNMAYG VKRLMSGPVP PSACPPRFSP ITIDGMTSLV GMNIPGHTGT GWCIFVYNLS PDSDESVLWQ LFGPFGAVNN VKVIRDFNTN KCKGFGFVTM TNYDEAAMAI ASLNGYRLGD RVLQVSFKTN KAHKS

image



Background

Alternative Names ELAV-like protein 4 (Hu-antigen D) (HuD) (Paraneoplastic encephalomyelitis)

Product Name: Recombinant Human ELAVL4 Protein
Catalog #: PHH8513



antigen HuD)

Note

For Research Use Only , Not for Diagnostic Use