


Product Name: Recombinant Human GNA12 Protein
Catalog #: PHH8540



Summary

Name	GNA12_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q03113
Host	HEK293
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	STFLKQMRIL HGREFDQKAL LEFRDTIFDN ILKGSRLVD ARDKLGIPWQ YSENEKHGMF LMAFENKAGL PVEPATFQLY VPALSALWRD SGIREAFSRR SEFQLGESVK YFLDNLDRIG QLNYFPSKQD ILLARKATKG IVEHDFVIKK

image



Background

Alternative Names Guanine nucleotide-binding protein subunit alpha-12 (G alpha-12) (G-protein subunit alpha-12)

Note

For Research Use Only , Not for Diagnostic Use