


Summary

Name	HERC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15751
Host	HEK293
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SRTCASPDSI QTGDAPIVSE TCEVYVWGSN SSHQLVEGTQ EKILQPKLAP SFSDAQTIEA GQYCTFVIST DGSVRACGKG SYGRLGLGDS NNQSTLKKLT FEPHRSIKKV SSSKGS DGHT LAF TTEGEVF SWGDGDY GKL GHGNSSTQKY

image



Background

Alternative Names	Probable E3 ubiquitin-protein ligase HERC1 (EC 2.3.2.26) (HECT domain and RCC1-like domain-containing protein 1) (HECT-type E3 ubiquitin transferase HERC1) (p532) (p619)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human HERC1 Protein
Catalog #: PHH8577

