


Product Name: Recombinant Human ARHGEF9 Protein
Catalog #: PHH8608



Summary

Name	ARHG9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	O43307
Host	HEK293
Species	Human
Molecular Mass	60 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SNKDWWWGQI DDEGWFPAS FVRLWVNQED EVEEGPSDVQ NGHLDPNSDC LCLGRPLQNR DQMRANVINE IMSTERHYIK HLKDICEGYL KQCRKRDMF SDEQLKVIFG NIEDIYRFQM GFVRDLEKQY NNDDPHLSEI GPCFLEHQDG

image



Background

Alternative Names	Rho guanine nucleotide exchange factor 9 (Collybistin) (PEM-2 homolog) (Rac/Cdc42 guanine nucleotide exchange factor 9)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use