Product Name: Recombinant Human ARHGEF9 Protein

Catalog #: PHH8609



Summary

Name ARHG9_HUMAN

Purity Greater than 85% as determined by reducing SDS-PAGE

Endotoxin level Please contact with the lab for this information

Tag C-His

Accession # 043307

Host HEK293

Species Human

Molecular Mass 19 kDa

Buffer Phosphate buffered saline

Shipping The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at

the temperature listed below.


Stability&Storage Store at \leq -70°C, stable for 6 months after receipt. Store at \leq -70°C, stable for 3

months under sterile conditions after opening. Please minimize freeze-thaw cycles. SNKDWWWGQI DDEEGWFPAS FVRLWVNQED EVEEGPSDVQ NGHLDPNSDC LCLGRPLQNR DQMRANVINE IMSTERHYIK HLKDICEGYL KQCRKRRDMF

SDEQLKVIFG NIEDIYRFQM GFVRDLEKQY NNDDPHLSEI GPCFLEHQDG

image

Sequence



Background

Alternative Names Rho guanine nucleotide exchange factor 9 (Collybistin) (PEM-2 homolog)

(Rac/Cdc42 guanine nucleotide exchange factor 9)

Note

For Research Use Only, Not for Diagnostic Use