


Product Name: Recombinant Human RBMS3 Protein
Catalog #: PHH8610




Summary

Name	RBMS3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-GST & C-His
Accession #	Q6XE24
Host	HEK293
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAKQQEQDPT NLYISNLPIS MDEQELENML KPFGHVISTR ILRDANGVSR GVGFARMEST EKCEVVIQHF NGKYLKTPPG IPAPSEPLLC KFADGGQKKR QNQSKYTQNG RPWPREGAG

image



SDS



WB

Background

Alternative Names RNA-binding motif, single-stranded-interacting protein 3

Note

For Research Use Only , Not for Diagnostic Use