



Product Name: Recombinant Human HERC6 Protein
Catalog #: PHH8905



Summary

Name	HERC6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q8IVU3
Host	HEK293
Species	Human
Molecular Mass	67 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DALCIRFSIH WDQNDVDLIP NGISIPVDQT NKRDYVSKYI DYIFNVSVKA VYEEFQRGFY RVCEKEILRH FYPEELMTAI IGNTDYDWKQ FEQNSKYEQG YQKSHPTIQL FWKAFHKLTL DEKKKFLFFL TGRDRLHARG IQKMEIVFRC PETFSERDHP TSITCHNILS LPKYSTMERM EEALQVAINN NRGFVSPMLT

image



Background

Alternative Names	Probable E3 ubiquitin-protein ligase HERC6 (EC 2.3.2.26) (HECT domain and RCC1-like domain-containing protein 6) (HECT-type E3 ubiquitin transferase HERC6)
--------------------------	---

Note

Product Name: Recombinant Human HERC6 Protein
Catalog #: PHH8905



For Research Use Only , Not for Diagnostic Use