



Product Name: Recombinant Human TDRD12 Protein
Catalog #: PHH8943



Summary

Name	TDR12_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q587J7
Host	HEK293
Species	Human
Molecular Mass	53 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LNAEMNEYFK DSNKTTVEKV EKFGLYGLAE KTLFHRVQVL EVNQKEDAWA LDDILVEFID EGRTGLVTRD QLLHLPEHFH

image



Background

Alternative Names Putative ATP-dependent RNA helicase TDRD12 (EC 3.6.4.13) (ES cell-associated transcript 8 protein) (Tudor domain-containing protein 12)

Note

For Research Use Only , Not for Diagnostic Use