



Summary

Name	UBP2L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q14157
Host	HEK293
Species	Human
Molecular Mass	58 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KKKGVSGQKD GGQTESNEEG KENRDRDRDY SRRGGPPRR GRGASRGREF RGQENGLDGT KSGGPSGRGT ERGRRGRGRG RGGSGRRGGR FSAQGMGTFN PADYAEPANT DDNYGNSSGN TWNNTGHFEP DDGTSAWRTA

image



Background

Alternative Names Ubiquitin-associated protein 2-like (Protein NICE-4)

Note

For Research Use Only , Not for Diagnostic Use