



Product Name: Recombinant Human WTIP Protein
Catalog #: PHH9075



Summary

Name	WTIP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	A6NIX2
Host	HEK293
Species	Human
Molecular Mass	59 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GRGRRGPRPG PGDEAAPALG RRGKGGSGGPE AGADGLSRGE RGPRAAVPE LSAQPAGSPR ASLAGSDGGG GGG SARSSGI SLGYDQRHGS PRSGRSDPRP GPGPPSVGSA RSSVSSLGSR GSAGAYADFL PPGACPAPAR SPEPAGPAPF PLPALPLPPG REGGPSAAER

image



Background

Alternative Names Wilms tumor protein 1-interacting protein (WT1-interacting protein)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human WTIP Protein
Catalog #: PHH9075

