



Product Name: Recombinant Human RND2 Protein
Catalog #: PHH9136



Summary

Name	RND2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P52198
Host	HEK293
Species	Human
Molecular Mass	67 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EGQSGRCKIV VVGDAECGKT ALLQVFAKDA YPGSYVPTVF ENYTASFEID KRRIELNMWD TSGSSYYDNV RPLAYPDSDA VLICFDISRP ETLDVLKKW QGETQEFCPN AKVVLVGCKL DMRTDLATLR ELSKQRLIPV THEQGTVLAK QVGAVSYVEC SSRSSERSVR DVFHVATVAS LGRGHRQLRR TDSRRGMQRS AQLSGRPDRG NEGEIHKDRA KSCNLM

image



Background

Alternative Names	Rho-related GTP-binding protein RhoN (Rho family GTPase 2) (Rho-related GTP-binding protein Rho7) (Rnd2)
--------------------------	--

Note

Product Name: Recombinant Human RND2 Protein
Catalog #: PHH9136



For Research Use Only , Not for Diagnostic Use