


Product Name: Recombinant Human RND2 Protein
Catalog #: PHH9137



Summary

Name	RND2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P52198
Host	HEK293
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NAKVVLVGCK LDMRTDLATL RELSKQRLIP VTHEQGTVLA KQVGAVSYVE CSSRSSERSV RDVFHVATVA SLGRGHRQLR RTDSRRGMQR SAQLSGRPDR GNEGEIHKDR

image



Background

Alternative Names	Rho-related GTP-binding protein RhoN (Rho family GTPase 2) (Rho-related GTP-binding protein Rho7) (Rnd2)
--------------------------	----------------------------------------------------------------------------------------------------------

Note

For Research Use Only , Not for Diagnostic Use