


Product Name: Recombinant Human MTMR14 Protein
Catalog #: PHH9300



Summary

Name	MTMRE_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NCE2
Host	HEK293
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ASAGSSASSG NQPPQELGLG ELLEEFSTRQ YRAKDGSGTG GSKVERIEKR CLELFGRDYC FSVIPNTNGD ICGHYPRHIV FLEYESSEKE KDTFESTVQV SKLQDLIHRS KMARCRGRFV CPVILFKGKH ICRSATLAGW GELYGRSGYN YFFSGGADDA WADVEDVTEE DCALRSGDTH LFDKVRGYDI

image



Background

Alternative Names Myotubularin-related protein 14 (EC 3.1.3.-) (HCV NS5A-transactivated protein 4 splice variant A-binding protein 1) (NS5ATP4ABP1) (hJumpy)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MTMR14 Protein
Catalog #: PHH9300

