



Product Name: Recombinant Human ARL14 Protein
Catalog #: PHH9507



Summary

Name	ARL14_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q8N4G2
Host	HEK293
Species	Human
Molecular Mass	64 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GSLGSKNPQT KQAQVLLGL DSAGKSTLLY KLKLAKDITT IPTIGFNVEM IELERNLSLT VWDVGGQEKM RTVWGCYCEN TDGLVYVVDS TDKQRLEESQ RQFEHILKNE HIKNVPVLL ANKQDMPGAL TAEDITRMFK VKKLCSDRNW YVQPCCALTG EGLAQGFRKL TGFVKSHMKS RGDTLAFFKQ N

image



Background

Alternative Names ADP-ribosylation factor-like protein 14 (ADP-ribosylation factor 7)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ARL14 Protein
Catalog #: PHH9507

