


Product Name: Recombinant Human BRAT1 Protein
Catalog #: PHH9535



Summary

Name	BRAT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q6PJG6
Host	HEK293
Species	Human
Molecular Mass	57 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QPVADDTCLE KLLDWFKTVT EGEVVLLQ EHPCLVELLS HVLKVQDLSS GVLSFSLRLA GTFAAQENCF QYLQQGELLP GLFGEPGPLG RATWAVPTVR SGWIIQLRSL AQHPSALRFL ADHGAVDTIF

image



Background

Alternative Names BRCA1-associated ATM activator 1 (BRCA1-associated protein required for ATM activation protein 1)

Note

For Research Use Only , Not for Diagnostic Use