


Product Name: Recombinant Human CELF4 Protein
Catalog #: PHH9578



Summary

| | |
|------------------------------|--|
| Name | CELF4_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | N-GST & C-His |
| Accession # | Q9BZC1 |
| Host | HEK293 |
| Species | Human |
| Molecular Mass | 40 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | IPRNLDEKDL KPLFEEFGKI YELTVLKDRF TGMHKGCAFL TYCERESALK AQSALHEQKT LPGMNRPIQV KPADESERGG SSCLRQPPSQ DRKLFVGMMLN KQQSEDDVRR |

image



Background

Alternative Names CUGBP Elav-like family member 4 (CELF-4) (Bruno-like protein 4) (CUG-BP- and ETR-3-like factor 4) (RNA-binding protein BRUNOL-4)

Note

For Research Use Only , Not for Diagnostic Use