


Product Name: Recombinant Human GAGE7 Protein
Catalog #: PHH9642



Summary

Name	GAGE7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	O76087
Host	HEK293
Species	Human
Molecular Mass	55 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSWRGRSTYY WPRPRRYVQP PEMIGPMRPE QFSDEVEPAT PEEGEPATQR QDPAAAEQEGE DEGASAGQGP KPEAHSQEQG HPQTGCECED GPDGQEMDPP NPEEVKTPEE GEKQSQC

image



Background

Alternative Names	G antigen 7 (GAGE-7) (AL4) (Cancer/testis antigen 4.7) (CT4.7) (GAGE-12I) (GAGE-7B) (GAGE-8)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use