


Product Name: Recombinant Human RIC8B Protein
Catalog #: PHH9916



Summary

Name	RIC8B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q9NVN3
Host	HEK293
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	THEETAQEAT TLDELPSNKT AEKETVLKNN TMVYNGMNME AIHVLLNFME KRIDKGSSYR EGLTPVLSLL TECSRAHRNI RKFLKDQVLP PLRDVTNRPE VGSTVRNKLKLV

image



Background

Alternative Names Synembryn-B (Brain synembryn) (hSyn) (Protein Ric-8B)

Note

For Research Use Only , Not for Diagnostic Use