


Summary

Name	Cathepsin D/CTSD
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Mouse Cathepsin D is produced by our Mammalian expression system and the target gene encoding Ile21-Leu410 is expressed with a 6His tag at the C-terminus.
Accession #	P18242
Host	Human Cells
Species	Mouse
Predicted Molecular Mass	43.9 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM MES, 150mM NaCl, pH 5.5.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. /xa0Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. /xa0

SDS-PAGE image

Product Name: Recombinant Mouse CSTD (C-6His)
Catalog #: PHM0243



Alternative Names

Cathepsin D; CTSD;CPSD

Background

CTSD localizes to the lysosome and consists of a light chain and a heavy chain. CTSD is expressed in epithelial cells as well as in macrophages. CTSD is a lysosomal aspartyl protease that depends critically on protonation of its active site Asp residue and gets activated at pH 5 in endosome of hepatocytes. It has been suggested to facilitate cancer cell migration and invasion by digesting the basement membrane, extracellular matrix and connective tissue. In addition, CTSD has been used as a breast cancer tumor marker.

Note

For Research Use Only , Not for Diagnostic Use.