

Product Name: Recombinant Mouse IL-9 (C-6His)
Catalog #: PHM1014




Summary

Name	IL-9/Interleukin-9
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Mouse Interleukin-9 is produced by our Mammalian expression system and the target gene encoding Gln19-Pro144 is expressed with a 6His tag at the C-terminus.
Accession #	P15247
Host	Human Cells
Species	Mouse
Predicted Molecular Mass	15.2 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Mouse IL-9 (C-6His)
Catalog #: PHM1014



Alternative Names

Interleukin-9; IL-9; Cytokine P40; T-Cell Growth Factor P40; IL9

Background

Interleukin-9 (IL-9) is a secreted protein that belongs to the IL-7/IL-9 family. Mature mouse IL-9 shares 57% and 74% amino acid sequence identity with human and rat IL-9, respectively. IL-9 supports IL-2 independent and IL-4 independent growth of helper T-cells. IL-9 stimulates cell proliferation and prevents apoptosis. It functions through the IL-9 receptor (IL-9R), which activates different signal transducer and activator (STAT) proteins and thus connects this cytokine to various biological processes. IL-9 has been identified as a candidate gene for asthma. IL-9 is a determining factor in the pathogenesis of bronchial hyperresponsiveness.

Note

For Research Use Only , Not for Diagnostic Use.