


Product Name: Recombinant Mouse Testican 3 (C-6His)
Catalog #: PHM1615



Summary

Name	Testican-3/SPOCK3
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Mouse Testican-3 is produced by our Mammalian expression system and the target gene encoding Ala23-Ile436 is expressed with a 6His tag at the C-terminus.
Accession #	Q8BKV0-1
Host	Human Cells
Species	Mouse
Predicted Molecular Mass	47.9 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image



Background

Product Name: Recombinant Mouse Testican 3 (C-6His)
Catalog #: PHM1615



Alternative Names

Testican-3;Spock3

Background

Testican3, also called SPOCK3, is a secreted protein and make up a family of extracellular heparan/chondroitin sulfate proteoglycans. It contains 1 Kazal-like domain and 1 thyroglobulin type-1 domain. Testican3 contain inhibitory regions in several domains targeted to different classes of protease, and in some cases may act as protease inhibitors. In addition to their presence in testis, testicans are enriched in brain and have been shown to regulate neuronal attachment and outgrowth. It may participate in diverse steps of neurogenesis and Inhibits the processing of pro-matrix metalloproteinase 2 (MMP-2) by MT1-MMP and MT3-MMP. It also may interfere with tumor invasion.

Note

For Research Use Only , Not for Diagnostic Use.