


Product Name: Recombinant Mouse APN Protein (C-6His)
Catalog #: PHM2232



Summary

Name	CD13/APN Protein/Aminopeptidase N
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Construction	Recombinant Mouse Aminopeptidase N is produced by our Mammalian expression system and the target gene encoding Lys69/xadSer966 is expressed with a 6His tag at the C-terminus.
Accession #	P97449
Host	Human Cells
Species	Mouse
Predicted Molecular Mass	103.6 KDa
Formulation	Supplied as a 0.2 µm filtered solution of PBS, pH 7.4.
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	

SDS-PAGE image



Background

Alternative Names	Aminopeptidase N; ANPEP; AP-M; APN; AP-N; CD13 antigen; CD13; CD13APN; PEPN; PEPNhAPN
Background	ANPEP gene encodes aminopeptidase N (APN) also known as microsomal aminopeptidase, alanyl aminopeptidase, aminopeptidase M, CD13, or membrane

Product Name: Recombinant Mouse APN Protein (C-6His)
Catalog #: PHM2232

protein p161, is a member of the peptidase M1 family. Widely expressed in many cells, tissues and species, APN cleaves the N-terminal amino acids from bioactive peptides, leading to their inactivation or degradation. Probably plays a role in regulating growth and differentiation of early B-lineage cells. It also may play a role in the catabolic pathway of the renin-angiotensin system. It degrades vasoconstricting angiotensin II into angiotensin III and therefore helps to regulate blood pressure.

Note

For Research Use Only , Not for Diagnostic Use.