

Product Name: Recombinant Vacciniavirus B18R (C-6His)
Catalog #: PHV1549




Summary

Name	Soluble Interferon α/β receptor B18
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/ μ g as determined by LAL test.
Construction	Recombinant Vacciniavirus Soluble Interferon Alpha/beta Receptor B18 is produced by our Mammalian expression system and the target gene encoding His20-Glu351 is expressed with a 6His tag at the C-terminus.
Accession #	P25213
Host	Human Cells
Species	Vaccinia virus
Predicted Molecular Mass	39.2 KDa
Formulation	Lyophilized from a 0.2 μ m filtered solution of PBS, pH 7.4.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Lyophilized protein should be stored at $\leq -20^{\circ}\text{C}$, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8 $^{\circ}\text{C}$ for 2-7 days. Aliquots of reconstituted samples are stable at $\leq -20^{\circ}\text{C}$ for 3 months.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μ g/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Vacciniavirus B18R (C-6His)
Catalog #: PHV1549



Alternative Names

Soluble interferon alpha/beta receptor B18; VACWR200

Background

B18R, also known as B19R in the Copenhagen strain of Vaccinia, is a type I interferon (IFN)-binding protein encoded by the B18R open reading frame in the Western Reserve strain of vaccinia virus. B18R exists in two forms, soluble and membrane-bound form and also has a broad species specificity. B18R has high affinity for human IFN-alpha and also binds rabbit, bovine, rat, pig, and mouse IFN-alpha and IFN-beta. Secreted B18R binds to uninfected and infected cells. B18R presents at the cell surface and protects cells from the antiviral state. It has shown that binding of soluble recombinant B18R protects cultured cells from IFN and allows vaccinia virus replication.

Note

For Research Use Only , Not for Diagnostic Use.