


Product Name: Recombinant Cynomolgus EDIL3 (C-6His)
Catalog #: PHV2438



Summary

Name	EDIL3/EGF-like repeats and discoidin I-like domains 3/Del-1/Integrin-binding protein DEL1
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Construction	Recombinant Cynomolgus EGF-like repeat and discoidin I-like domain-containing protein 3 is produced by our Mammalian expression system and the target gene encoding Asp24-Glu480 is expressed with a 6His tag at the C-terminus.
Accession #	A0A2K5WYU3
Host	Human Cells
Species	Cynomolgus
Predicted Molecular Mass	52.1 kDa
Formulation	Lyophilized protein should be stored at $\leq -20^{\circ}\text{C}$, stable for one year after receipt. Reconstituted protein solution can be stored at $2-8^{\circ}\text{C}$ for 2-7 days. Aliquots of reconstituted samples are stable at $\leq -20^{\circ}\text{C}$ for 3 months.
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	0.0
Reconstitution	

SDS-PAGE image



Background

Alternative Names	0.0
--------------------------	-----

Product Name: Recombinant Cynomolgus EDIL3 (C-6His)
Catalog #: PHV2438



Background

0.0

Note

For Research Use Only , Not for Diagnostic Use.