

Production Name: R-PE-conjugated Monoclonal Mouse Anti-Chicken IgM, μ Chain Specific Secondary Antibody

Catalog #: AMS0498

Note:Centrifuge before opening to ensure complete recovery of vial contents

Summary

| | |
|------------------------|--|
| Production Name | R-PE-conjugated Monoclonal Mouse Anti-Chicken IgM, μ Chain Specific Secondary Antibody |
| Clonality | Monoclonal IgG1,Kappa |
| Host | Mouse |
| Reactivity | Chicken |
| Application | IF,FC |
| Form | Liquid |
| Concentration | 0.5mg/mL |
| Specificity | IgM, μ Chain Specific |
| Buffer | 0.01M PBS,0.5%protective protein,0.1%Proclin-300,pH7.4 |
| Conjugation | R-PE |
| Ex-max | 496nm,565nm |
| Em-max | 575nm |
| Purification | Affinity purification |
| Purity | Calculated as antibodies, no less than 90%. |
| Storage | Store at 2-8°C and do not freeze. Prepare working dilution on day of use.Stable for more than six months after shipment. |

Application

1:100-1:1000.This data is a recommended range, please refer to the actual value according to the specific situation. The optimal dilution ratio is influenced by various factors, such as antigen concentration and permeability. The actual dilution ratio to be used must be determined through experiments.

Note

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.