

**Production Name: Monoclonal Mouse Anti-Monkey IgG Secondary Antibody****Catalog #: AMS0557**

---

Note:Centrifuge before opening to ensure complete recovery of vial contents

**Summary**

<b>Production Name</b>	Monoclonal Mouse Anti-Monkey IgG Secondary Antibody
<b>Clonality</b>	Monoclonal IgG1,Kappa
<b>Host</b>	Mouse
<b>Reactivity</b>	Monkey
<b>Application</b>	ELISA,IF,FC,WB,IHC
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/mL
<b>Specificity</b>	IgG
<b>Buffer</b>	0.01M PBS,pH7.2,0.02% N preservative
<b>Conjugation</b>	Unconjugated
<b>Ex-max</b>	none
<b>Em-max</b>	none
<b>Purification</b>	Affinity purification
<b>Purity</b>	Calculated as antibodies, no less than 90%.
<b>Storage</b>	Store at -20°C.Avoid repeated freezing and thawing. Stable for more than one year after shipment.

**Application**

Used in WB,ELISA,IF,IP,IHC,FCM,etc.Before labeling, perform dialysis or ultrafiltration to replace the solution.

**Note**

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.