

Production Name: HRP-conjugated Monoclonal Mouse Anti-Dog IgM, μ Chain Specific Secondary Antibody

Catalog #: AMS0573

Note:Centrifuge before opening to ensure complete recovery of vial contents

Summary

Production Name	HRP-conjugated Monoclonal Mouse Anti-Dog IgM, μ Chain Specific Secondary Antibody
Clonality	Monoclonal IgG2b,Kappa
Host	Mouse
Reactivity	Dog
Application	ELISA, WB, Dot blot
Form	Liquid
Concentration	0.5mg/mL
Specificity	IgM, μ Chain Specific
Buffer	50% Glycerol,0.01M PBS,0.5%protective protein,0.05%Proclin-300,pH7.4
Conjugation	HRP
Ex-max	none
Em-max	none
Purification	Affinity purification
Purity	Calculated as antibodies, no less than 90%.
Storage	Please store at -20°C away from light.Stable for one year after shipment.

Application

WB 1:2000-10000(ECL substrates)

WB 1:1000-5000(DAB substrates)

Direct ELISA 1:5000-30000(TMB substrates)

IHC 1:500-2000(DAB substrates).This data is a recommended range, please refer to the actual value according to the specific situation. The optimal dilution ratio is influenced by various factors, such as antigen concentration and permeability. The actual dilution ratio to be used must be determined through experiments.

Note

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.