

Production Name: R-PE-conjugated Polyclonal Goat Anti-Rabbit IgG(H+L)(Non cross with Mouse IgG) Secondary Antibody

Catalog #: APS0155

Note:Centrifuge before opening to ensure complete recovery of vial contents

Summary

Production Name	R-PE-conjugated Polyclonal Goat Anti-Rabbit IgG(H+L)(Non cross with Mouse IgG) Secondary Antibody
Clonality	Polyclonal Whole IgG
Host	Goat
Reactivity	Rabbit
Application	IF,FC
Form	Liquid
Concentration	0.5mg/mL
Specificity	IgG(H+L)(Non cross with Mouse IgG)
Buffer	0.01M PBS,0.5%protective protein,0.1%Proclin-300,pH7.4
Conjugation	R-PE
Ex-max	496nm,565nm
Em-max	575nm
Purification	Affinity purification
Purity	Calculated as antibodies, no less than 90%.
Storage	Store at 2-8°C and do not freeze. Prepare working dilution on day of use.Stable for more than six months after shipment.

Application

1:100-1:1000.This data is a recommended range, please refer to the actual value according to the specific situation. The optimal dilution ratio is influenced by various factors, such as antigen concentration and permeability. The actual dilution ratio to be used must be determined through experiments.

Note

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.