

**Production Name:** Polyclonal Goat Anti-Mouse IgG(H+L)(Non cross with Rabbit IgG)  
**Secondary Antibody**

**Catalog #:** APS0641

---

Note:Centrifuge before opening to ensure complete recovery of vial contents

## Summary

|                        |                                                                                                   |
|------------------------|---------------------------------------------------------------------------------------------------|
| <b>Production Name</b> | Polyclonal Goat Anti-Mouse IgG(H+L)(Non cross with Rabbit IgG) Secondary Antibody                 |
| <b>Clonality</b>       | Polyclonal Whole IgG                                                                              |
| <b>Host</b>            | Goat                                                                                              |
| <b>Reactivity</b>      | Mouse                                                                                             |
| <b>Application</b>     | ELISA,IF,FC,WB,IHC                                                                                |
| <b>Form</b>            | Liquid                                                                                            |
| <b>Concentration</b>   | 1mg/mL                                                                                            |
| <b>Specificity</b>     | IgG(H+L)(Non cross with Rabbit IgG)                                                               |
| <b>Buffer</b>          | 0.01M PBS,pH7.2,0.02% N preservative                                                              |
| <b>Conjugation</b>     | Unconjugated                                                                                      |
| <b>Ex-max</b>          | none                                                                                              |
| <b>Em-max</b>          | none                                                                                              |
| <b>Purification</b>    | Affinity purification                                                                             |
| <b>Purity</b>          | Calculated as antibodies, no less than 90%.                                                       |
| <b>Storage</b>         | Store at -20°C.Avoid repeated freezing and thawing. Stable for more than one year after shipment. |

## Application

Used in WB,ELISA,IF,IP,IHC,FCM,etc.Before labeling, perform dialysis or ultrafiltration to replace the solution.

## Note

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.