

**Production Name:** Polyclonal Goat Anti-Rabbit IgG(H+L)(Non cross with Mouse IgG)  
**Secondary Antibody**

**Catalog #:** APS0645

---

Note:Centrifuge before opening to ensure complete recovery of vial contents

## Summary

<b>Production Name</b>	Polyclonal Goat Anti-Rabbit IgG(H+L)(Non cross with Mouse IgG) Secondary Antibody
<b>Clonality</b>	Polyclonal Whole IgG
<b>Host</b>	Goat
<b>Reactivity</b>	Rabbit
<b>Application</b>	ELISA,IF,FC,WB,IHC
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/mL
<b>Specificity</b>	IgG(H+L)(Non cross with Mouse IgG)
<b>Buffer</b>	0.01M PBS,pH7.2,0.02% Sodium Azide
<b>Conjugation</b>	Unconjugated
<b>Ex-max</b>	none
<b>Em-max</b>	none
<b>Purification</b>	Affinity purification
<b>Purity</b>	Calculated as antibodies, no less than 90%.
<b>Storage</b>	Store at -20°C.Avoid repeated freezing and thawing. Stable for more than one year after shipment.

## Application

Used in WB,ELISA,IF,IP,IHC,FCM,etc.Before labeling, perform dialysis or ultrafiltration to replace the solution.

## Note

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.