

Production Name: Polyclonal Goat Anti-Mouse IgG(H+L)(Non cross with Human IgG)
Secondary Antibody

Catalog #: APS0649

Note:Centrifuge before opening to ensure complete recovery of vial contents

Summary

Production Name	Polyclonal Goat Anti-Mouse IgG(H+L)(Non cross with Human IgG) Secondary Antibody
Clonality	Polyclonal Whole IgG
Host	Goat
Reactivity	Mouse
Application	ELISA,IF,FC,WB,IHC
Form	Liquid
Concentration	1mg/mL
Specificity	IgG(H+L)(Non cross with Human IgG)
Buffer	0.01M PBS,pH7.2,0.02% N preservative
Conjugation	Unconjugated
Ex-max	none
Em-max	none
Purification	Affinity purification
Purity	Calculated as antibodies, no less than 90%.
Storage	Store at -20°C.Avoid repeated freezing and thawing. Stable for more than one year after shipment.

Application

Used in WB,ELISA,IF,IP,IHC,FCM,etc.Before labeling, perform dialysis or ultrafiltration to replace the solution.

Note

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.