

Production Name: Alkaline Phosphatase-conjugated Streptavidin

Catalog #: RAD25002

Note:Centrifuge before opening to ensure complete recovery of vial contents

Summary

Production Name	Alkaline Phosphatase-conjugated Streptavidin
Reactivity	Biotin
Application	IHC/ICC,ELISA,WB
Form	Liquid
Concentration	0.5mg/mL
Buffer	50% Glycerol,0.005 M Tris-HCl, 0.125 M NaCl, 0.5%protective protein,0.02% N preservative,pH 8.0
Conjugation	Alkaline Phosphatase
Ex-max	none
Em-max	none
Storage	Please store at -20°C and avoid repeated freezing and thawing .Stable for one year after shipment.

Application

IHC/ICC:1:250-500

ELISA、WB:1:500-5000.This data is a recommended range, please refer to the actual value according to the specific situation. The optimal dilution ratio is influenced by various factors, such as antigen concentration and permeability. The actual dilution ratio to be used must be determined through experiments.

Note

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this product into small packages and use one tube each time. This product is only for research use, not for diagnosis.