

Production Name: R-PE-conjugated Streptavidin

Catalog #: RAD25015

Note:Centrifuge before opening to ensure complete recovery of vial contents

Summary

Production Name	R-PE-conjugated Streptavidin
Reactivity	Biotin
Application	FC,Immunoassay
Form	Liquid
Concentration	0.5mg/mL
Buffer	0.01M PBS,0.5%protective protein,0.1%Proclin-300,pH 7.4
Conjugation	R-PE
Ex-max	496nm,565nm
Em-max	575nm
Storage	Please store at 2-8°C away from light and absolutely avoid freezing.Stable for more than 6 months after shipment.

Application

0.1-2 µg/mL (For most applications, Most commonly used in flow cytometry applications) .This data is a recommended range, please refer to the actual value according to the specific situation. The optimal dilution ratio is influenced by various factors, such as antigen concentration and permeability. The actual dilution ratio to be used must be determined through experiments.

Note

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this product into small packages and use one tube each time. This product is only for research use, not for diagnosis.