

Production Name:Flag-Tag Mouse Monoclonal Antibody
Catalog#: AMM80002

Summary

Production Name	Flag-Tag Mouse Monoclonal Antibody
Host	Mouse
Application	WB,Elisa,IP,IF,IHC,ICC
specificity	The antibody detects C-terminal, internal, and N-terminal Flag-tag fusion proteins.
Reactivity	Species independent
Isotype	IgG1
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	PBS,PH:7.2, 50% Glycerol
Purification	Affinity purification by Protein A

Immunogen Information

Gene Name
Immunogen

Application

Dilution Ratio	WB 1:500-1:5000, ELISA 1:5000-20000, IP 1:200-500, IF-P/IF-F/IF-ICC 1:200-500, IHC 1:100-400
-----------------------	--

Note

For research use only .