

Production Name:GFP-Tag Monoclonal Antibody(1G6)
Catalog#: AMM80005



Summary

Production Name	GFP-Tag Monoclonal Antibody(1G6)
Host	Mouse
Application	WB,Elisa
specificity	
Reactivity	Species independent
Isotype	IgG1
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 2% BSA and New type preservative N.
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.

Immunogen Information


Gene Name
Immunogen

Application

Dilution Ratio WB 1:5000, ELISA 1:5000-20000, Optimal dilutions should be determined by the end user

Background

Image Data



2ug GFP fusion protein+ Primary antibody dilution at 1) 1:5000;2) 1:10000

Production Name:GFP-Tag Monoclonal Antibody(1G6)
Catalog#: AMM80005



Note

For research use only .