

Product Name: HA-Tag Monoclonal Antibody(1B10)**Catalog #: AMM80008**

For research use only.

Summary


Host	Mouse
Application	WB,IP,IF
Specificity	The antibody detects HA-tag fusion proteins.
Reactivity	Species independent
Isotype	IgG
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 2% BSA and New type preservative N.
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.

Application

Dilution Ratio	WB 1:5000, IP 1:200, ICC/IF 1:1000, ELISA 1:5000-20000
-----------------------	--

Antigen Information

Gene Name	HA-tag
Immunogen	

Background**Image Data**

0.5ug HA fusion protein+ Primary antibody dilution at 1:10000

