

Production Name:HA-Tag Monoclonal Antibody(1B10)
Catalog#: AMM80008



Summary

Production Name	HA-Tag Monoclonal Antibody(1B10)
Host	Mouse
Application	WB,IP,IF
specificity	
Reactivity	Species independent
Isotype	
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 2% BSA and New type preservative N.
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.

Immunogen Information


Gene Name
Immunogen

Application

Dilution Ratio WB 1:5000, IP 1:200, IF-P/IF-F/IF-ICC 1:1000, ELISA 1:5000-20000

Background

Image Data



0.5ug HA fusion protein+ Primary antibody dilution at 1:10000

Production Name:HA-Tag Monoclonal Antibody(1B10)
Catalog#: AMM80008



Note

For research use only .