

Product Name: His-Tag Monoclonal Antibody(4E6)**Catalog #: AMM80010**

For research use only.

Summary

Host	Mouse
Application	WB,IP,IF
Specificity	The antibody detects C-terminal, internal, and N-terminal His-tag fusion proteins.
Reactivity	Species independent
Isotype	IgG1
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 2% BSA and New type preservative N.
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.


Application

Dilution Ratio	WB 1:3000, IP 1:200, ICC/IF 1:1000, ELISA 1:5000-20000
-----------------------	--

Antigen Information

Gene Name	His-tag
Immunogen	Synthetic peptide:HHHHHH conjugated to KLH.

Background**Image Data**



2ug His fusion protein+ Primary antibody dilution at 1) 1:5000,2) 1:10000