

Production Name: His-Tag Monoclonal Antibody(4E6)
Catalog#: AMM80010



Summary

Production Name	His-Tag Monoclonal Antibody(4E6)
Host	Mouse
Application	WB,IP,IF
specificity	Synthetic peptide:HHHHHH conjugated to KLH.
Reactivity	Species independent
Isotype	IgG1
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 2% BSA and New type preservative N.
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.

Immunogen Information

Gene Name
Immunogen


Application

Dilution Ratio WB 1:3000, IP 1:200, IF-P/IF-F/IF-ICC 1:1000, ELISA 1:5000-20000

Background

Image Data

Production Name:His-Tag Monoclonal Antibody(4E6)
Catalog#: AMM80010



2ug His fusion protein+ Primary antibody dilution at 1) 1:5000;2) 1:10000

Note

For research use only .