

Production Name: Plant actin Monoclonal Antibody(Q30)
Catalog#: AMM80015

Summary

Production Name	Plant actin Monoclonal Antibody(Q30)
Host	Mouse
Application	WB, Elisa
specificity	The antibody detects endogenous plant actin protein.
Reactivity	Various Plants
Isotype	
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS, pH 7.4, containing 50% glycerol, 2% BSA and New type preservative N.
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.

Immunogen Information

Gene Name	
Immunogen	Recombinant Protein of Plant actin

Application


Dilution Ratio	WB 1:2000-5000
-----------------------	----------------

Background

Actin is a highly conserved protein and an essential component of cell cytoskeleton and plays an important role in cytoplasmic streaming, cell shape determination, cell division, organelle movement and extension growth

Image Data

Production Name: Plant actin Monoclonal Antibody(Q30)
Catalog#: AMM80015



Western blot analysis of Arabidopsis, (Q30) diluted at 1:5000

Note

For research use only .