

Product Name: V5-Tag Monoclonal Antibody(3C8)**Catalog #: AMM80016**

For research use only.

Summary

Host	Mouse
Application	WB,IP,IF
Specificity	The antibody detects C-terminal, internal, and N-terminal V5-tag fusion proteins.
Reactivity	Species independent
Isotype	IgG1
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 2% BSA and New type preservative N.
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.

Application

Dilution Ratio	WB 1:5000, IP 1:200, IF-P/IF-F/IF-ICC 1:1000, ELISA 1:5000-20000
-----------------------	--

Antigen Information

Gene Name	V5-Tag
Immunogen	Synthetic Peptide of V5-Tag

Background**Image Data**