

Production Name:V5-Tag Monoclonal Antibody(3C8)
Catalog#: AMM80016



Summary

Production Name	V5-Tag Monoclonal Antibody(3C8)
Host	Mouse
Application	WB,IP,IF
specificity	The antibody detects C-terminal, internal, and N-terminal V5-tag fusion proteins.
Reactivity	Species independent
Isotype	IgG1
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 2% BSA and New type preservative N.
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.

Immunogen Information


Gene Name	
Immunogen	Synthetic Peptide of V5-Tag

Application

Dilution Ratio	WB 1:5000, IP 1:200, IF-P/IF-F/IF-ICC 1:1000, ELISA 1:5000-20000
-----------------------	--

Background

Image Data



1ug V5 fusion protein+ Primary antibody dilution at 1) 1:5000; 2) 1:10000

Production Name:V5-Tag Monoclonal Antibody(3C8)
Catalog#: AMM80016



Note

For research use only .