

Production Name:β-actin Monoclonal Antibody(5B7)
Catalog#: AMM80019



Summary

Production Name	β-actin Monoclonal Antibody(5B7)
Host	Mouse
Application	IF, WB, IHC-p, Elisa
specificity	
Reactivity	Human, Rat, Mouse, Mk, Dg, Ch, Hamster, Rb, Insect
Isotype	
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS, pH 7.4, containing 50% glycerol, 2% BSA and New type preservative N.
Purification	The antibody detects endogenous β-actin protein.

Immunogen Information


Gene Name	
Immunogen	Synthetic Peptide of β-actin

Application

Dilution Ratio	IF-P/IF-F/IF-ICC 1:50-200, WB 1:5000, IHC-P 1:200
-----------------------	---

Background

Image Data



Western blot analysis of 1 HEK293 2 SW480 3 HEPG2 4 MCF-7 5 mouse brain 6 Rat brain 7 HeLa 8 A549 lysates, primary antibody was diluted at 1:5000, 4° over night, secondary antibody(cat: RS23910) was diluted at 1:10000, 37° 1 hour.

Production Name:β-actin Monoclonal Antibody(5B7)
Catalog#: AMM80019



Note

For research use only .