
Product Name: Trx Tag Mouse Monoclonal Antibody**Catalog #:** AMM80032

For research use only.

Summary

Host	Mouse
Application	WB
Specificity	The antibody detects Trx-tag fusion proteins.
Reactivity	Species-independent
Isotype	IgG1
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB-1:1000-1:10000
-----------------------	-------------------


Antigen Information

Gene Name	Trx
Immunogen	Recombinant Protein of Trx-Tag

Background

Thioredoxin (Trx) is a redox protein that is found in several species, such as bacteria, plants and mammals, and contains a conserved active site, consisting of Trp-Cys-Gly-Pro-Cys. Participates in many cellular processes including redox signaling, response to oxidative stress, and protein reduction.

Image Data



1ug Trx fusion protein+ Primary antibody dilution at 1) 1:5000. 2) 1:10000