

Product Name: E2 Tag Mouse Monoclonal Antibody**Catalog #: AMM80042**

For research use only.

Summary

Host	Mouse
Application	WB
Specificity	The antibody detects E2-tag fusion proteins.
Reactivity	Species-independent
Isotype	IgG1
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB-1:1000-1:10000
-----------------------	-------------------


Antigen Information

Gene Name	E2
Immunogen	Synthetic Peptide of E2-Tag

Background

The E2 Tag is a peptide which is derived from the Bovine Papillomavirus type-1 transactivator protein E2 and consists of 10 amino acids (SSTSSDFRDR). The detection of E2 tagged proteins is based on the binding of mouse monoclonal antibodies specific to the tagged recombinant protein.

Image Data



2ug E2 fusion protein+ Primary antibody dilution at 1) 1:5000. 2) 1:10000