
Product Name: S Tag Mouse Monoclonal Antibody**Catalog #: AMM80043**

For research use only.

Summary

Host	Mouse
Application	WB
Specificity	The antibody detects S-tag fusion proteins.
Reactivity	Species-independent
Isotype	IgG1
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB-1:1000-1:10000
-----------------------	-------------------


Antigen Information

Gene Name	S
Immunogen	Synthesized peptide of S-tag

Background

S tag is an epitope tag composed of a 15 residue peptide, KETAAAKFERQHMD_S, derived from the pancreatic ribonuclease A.

Image Data



Western blot analysis of Stag (8G4) in Multitag antigen 120ng, Multitag antigen 210ng lysates using Stag (8G4) antibody.